ROTO SILENTA 630 RS

Ⓧ DE	Bedienungsanleitung.....	5
Ⓧ EN	Operating Instructions.....	27
Ⓧ FR	Mode d'emploi	48
Ⓧ IT	Istruzioni per l'uso	70


Fig. 1


Fig. 2

CE 0482

**EG-Konformitätserklärung
EC Conformity Declaration
Déclaration de conformité CE
Dichiarazione di conformità alle norme CEE**

Andreas Hettich GmbH & Co. KG • Föhrenstraße 12 • D-78532 Tuttlingen • Germany

Das bezeichnete Gerät, inklusive Zubehör entspricht den aufgeführten EG-Richtlinien und Normen.
The denoted device, including accessories corresponds to the listed EC guidelines and standards.
L'appareil désigné, y compris les accessoires, correspond aux directives CE et aux normes énumérées.
L'apparecchio designato, compresi gli accessori, è conforme alle direttive CE e alle norme citate.

Geräteart, Type of device, Type d'appareil, Tipo di apparecchio:

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Contents

1	Use according to specification	29
2	Residual risks	29
3	Technical specifications	29
4	Notes on safety	31
5	Symbol meanings	32
6	Delivery checklist	32
7	Unpacking the centrifuge and installation	33
8	Initial operation	33
9	Opening and closing the lid	33
9.1	Opening the lid	33
9.2	Closing the lid	33
10	Installation and removal of the rotor	34
11	Loading the rotor	34
12	Control and display devices	34
12.1	Status symbols	34
12.2	Control knob	35
12.3	Keys and set-up options	35
13	Entering centrifugation parameters	36
13.1	Running time	36
13.1.1	Time selection	36
13.1.2	Continuous running	36
13.2	Speed (RPM)	36
13.3	Relative centrifugal force (RCF/RZB)	36
13.4	Run-up and run-down parameters	36
13.4.1	Run-up step	36
13.4.2	Run-up time	36
13.4.3	Braking step	36
13.4.4	Run-down time	36
13.4.5	Brake de-energisation speed	37
13.5	Radius/Temperature	37
13.5.1	Radius	37
13.5.2	Temperature (for centrifuges with cooling only)	37
13.6	Automatic buffer	37
14	Programming	37
14.1	Entering / modifying programs	37
14.2	Calling up programs	37
15	Centrifugation	37
15.1	Centrifugation with time selection	38
15.2	Continuous running	38
16	Changing settings during the centrifugation run	38
17	Integral RCF (IRCF)	38
18	Displaying the maximum speed of the rotor	38

19	Displaying the maximum RCF of the rotor.....	38
20	Emergency stop	38
21	Audible signal.....	39
22	Interrogating the operating hours	39
23	Setting the date and time	39
24	Centrifugation data displayed after switching on.....	39
25	Immediate display of the centrifugation data after switching on	39
26	Key-operated switch (for centrifuges with key-operated switch only).....	40
27	Program linking (for centrifuges with program linking only).....	40
27.1	Linking programs	40
27.2	Changing program links	40
27.3	Centrifugation run with program linking.....	40
27.4	Deleting program linking	41
28	Cooling (only in centrifuges with cooling).....	41
28.1	Standby-cooling	41
28.2	Pre-cooling the rotor	41
29	Heating (only on centrifuges with heating/cooling option)	41
30	Relative centrifugal force (RCF).....	41
31	Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm ³	42
32	Rotor Identification	42
33	Emergency release	42
34	Maintenance and servicing.....	43
34.1	Centrifuge (housing, lid and centrifuging chamber).....	43
34.1.1	Surface cleaning and care	43
34.1.2	Surface disinfection	43
34.1.3	Removal of radioactive contaminants	43
34.2	Rotors and Attachments	44
34.2.1	Cleaning and care.....	44
34.2.2	Disinfection	44
34.2.3	Removal of radioactive contaminants	44
34.2.4	Trunnions.....	44
34.2.5	Rotors and accessories with limited service lives	44
34.3	Autoclaving	45
34.4	Centrifuge containers.....	45
35	Faults	46
36	Acceptance of the centrifuges for repair.....	47
37	Disposal	47
38	Anhang / Appendix.....	92
38.1	Rotoren und Zubehör / Rotors and accessories.....	92

1 Use according to specification

This machine is a medical device (laboratory centrifuge) in the sense of IVD guideline 98/79/EC and the guideline for medical products 93/42/EEC. The centrifuge is used to separate substances or substance mixtures with a density of max. 1.2 kg/dm³. This also includes substances and substance mixtures of human origin. The centrifuge is only intended to be used for this purpose. A different use or application over and above this is deemed not in accordance with the specifications. The company Andreas Hettich GmbH & Co. KG undertakes no liability for damages resulting therefrom.

Belonging to the application according to specification is also the observance of all references contained in the Instruction Manual and compliance with the inspection and maintenance works.

2 Residual risks

The machine is constructed according to the state of the art and the recognized technical safety regulations. Improper use and handling can result in dangers to life and limb of the user or third parties and impairments to the machine or to other material assets. The machine is only to be used for the specified applications and only in an impeccable technical safety condition.

Disturbances that can interfere with the safety are to be immediately rectified.

3 Technical specifications

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen			
Model	ROTO SILENTA 630 RS			
Type	5005, 5005-20, 5005-50, SK 48.06-30	5005-80	5005-90	5005-60
Mains voltage ($\pm 10\%$)	400 V 3~ + N			
Mains frequency	50 - 60 Hz			
Connected load	9700 VA	9400 VA	6600 VA	
Current consumption	14 A	13.5 A	9.5 A	
Cooling medium	R 404A			----
Max. capacity	12000 ml			
Allowed density	1.2 kg/dm ³			
Speed (RPM)	4500			
Force (RCF)	6498			
Kinetic energy	215000 Nm			
Obligatory inspection (BGR 500)	Yes			
Ambient conditions (EN / IEC 61010-1)	Indoors only Up to 2000 m above sea level 5°C to 40°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.			
– Set-up site				
– Altitude				
– Ambient temperature				
– Humidity				
– Excess-voltage category (IEC 60364-4-443)	II			
– Pollution degree	2			
Device protection class	I			
Not suitable for use in explosion-endangered areas.				
EMC	EN / IEC 61326-1, Class B			
– Emitted interference, Interference immunity				
Noise level (dependent on rotor)	≤ 62 dB(A)		≤ 56 dB(A)	
Dimensions	813 mm			
– Width				
– Depth	1015 mm	1050 mm		
– Height	973 mm			
Weight	approx. 360 kg	approx. 367 kg	approx. 306 kg	approx. 303 kg

Manufacturer	Andreas Hettich GmbH & Co. KG D-78532 Tuttlingen			
Model	ROTO SILENTA 630 RS			
Type	5005-08A	5005-08	5005-88	5005-98
Mains voltage ($\pm 10\%$)	220 V +6/-10% 3~ (+N)+PE	208 V 3~ + N		
Mains frequency	60 Hz			
Connected load	9000 VA		8700 VA	5200 VA
Current consumption	25 A		24 A	14 A
Cooling medium	R 404A			
Max. capacity	12000 ml			
Allowed density	1.2 kg/dm ³			
Speed (RPM)	4500			
Force (RCF)	6498			
Kinetic energy	215000 Nm			
Obligatory inspection (BGR 500)	yes			
Ambient conditions (EN / IEC 61010-1)	<p>Indoors only Up to 2000 m above sea level</p> <p>5°C to 30°C 5°C to 40°C Maximum relative humidity 80% for temperatures up to 31°C, linearly decreasing to 50% relative humidity at 40°C.</p>			
– Set-up site				
– Altitude				
– Ambient temperature				
– Humidity				
– Excess-voltage category (IEC 60364-4-443)	II			
– Pollution degree	2			
Device protection class	I			
Not suitable for use in explosion-endangered areas.				
EMC	EN / IEC 61326-1, Class B			
– Emitted interference, Interference immunity				
Noise level (dependent on rotor)	≤ 62 dB(A)			
Dimensions				
– Width	813 mm			
– Depth	1015 mm			
– Height	973 mm			
Weight	approx. 401 kg			

4 Notes on safety

No claim of warranty will be considered by the manufacturer unless ALL instructions in this manual have been followed.

- **The centrifuge should be installed on a good, stable base.**
- **Before using the centrifuge absolutely check the rotor for firm placement.**
- **When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.**
- **Rotors, suspensions and accessories that possess traces of corrosion or mechanical damage or if their term of use has expired may not be used any longer.**
- **The centrifuge may no longer be put into operation when the centrifuging chamber has safety-related damages.**
- **With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.**

- **Before the initial operation of your centrifuge you should read and pay attention to the operating instructions. Only personnel that has read and understood the operating instructions are allowed to operate the device.**
- Along with the operating instructions and the legal regulations on accident prevention, you should also follow the recognised professional regulations for working in a safe and professional manner. These operating instructions should be read in conjunction with any other instructions concerning accident prevention and environmental protection based on the national regulations of the country where the device is to be used.
- This centrifuge is a state-of-the-art piece of equipment which is extremely safe to operate. However, it can lead to danger for users or others if used by untrained staff, in an inappropriate way or for a purpose other than that it was designed for.
- The centrifuge must not be moved or knocked during operation.
- In case of fault or emergency release, never touch the rotor before it has stopped turning.
- To avoid damage due to condensate, when changing from a cold to a warm room the centrifuge must either heat up for at least 3 hours in the warm room before being connected to the mains, or run hot for 30 minutes in the cold room.
- Only the rotors and accessories approved by the manufacturer for this device may be used (see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").
- The centrifuge rotor may only be loaded in accordance with the chapter "Loading the rotor".
- When centrifuging with maxim revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.
- The centrifuge may only be operated when the balance is within the bounds of acceptability.
- The centrifuge may not be operated in explosion-endangered areas.
- The centrifuge must not be used with:
 - inflammable or explosive materials
 - materials that react with one another producing a lot of energy.
- If users have to centrifuge hazardous materials or compounds contaminated with toxic, radioactive or pathogenic micro-organisms, they must take appropriate measures.
 For hazardous substances centrifuge containers with special screw caps must strictly be used. In addition to the screw cap centrifuge containers, for materials in hazard category 3 and 4 a biosafety system must be used (see the World Health Organisation's "Laboratory Biosafety Manual").
 In a biosafety system, droplets and aerosols are prevented from escaping by a bioseal (packing ring).
 If the hanger of a biosafety system is used without the lid, the packing ring must be removed from the hanger in order to prevent the packing ring from being damaged during the centrifugation run. Damaged packing rings must not be used to seal the biosafety system.
 Without the use of a biosafety system the centrifuge is not microbiologically sealed in the sense of the EN / IEC 610101-2-020 standard.

For further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". If in doubt, you should obtain relevant information from the manufacturer.

- The centrifuge must not be operated with highly corrosive substances which could impair the mechanical integrity of rotors, hangers and accessories.
- Repairs must only be carried out by personnel authorised to do so by the manufacturer.
- Only original spare parts and original accessories licensed by the Andreas Hettich GmbH & Co. KG company are allowed to be utilised.
- The following safety regulations apply:
EN / IEC 61010-1 and EN / IEC 61010-2-020 as well as their national deviations.
- The safe operation and reliability of the centrifuge can only be guaranteed if:
 - the centrifuge is operated in accordance with the operating instructions,
 - the electrical installation on the site where the centrifuge is installed conforms to the demands of EN / IEC stipulations,
 - prescribed tests to BGV A1, BGR 500 are carried out by an expert.

5 Symbol meanings

Symbol on the machine:
Attention, general hazard area.
Before using the centrifuge implicitly read the operating instructions and pay attention to the safety relevant references!

Symbol in the operating instructions:
Attention, general hazard area.
This symbol refers to safety relevant warnings and indicates possibly dangerous situations.
The non-adherence to these warnings can lead to material damage and injury to personal.

Symbol on the machine and in the operating instructions:
Beware of hot surface.

Symbol on the machine:
Equipotential: Connector (PE connector) for potential equalization (only for centrifuge with PE connector).

Symbol on the machine:
Switch positions of the key-operated switch.

Symbol in the operating instructions:
This symbol refers to important circumstances.

Symbol on the machine and in the operating instructions:
Symbol for the separate collection of electric and electronic devices according to the guideline 2002/96/EG (WEEE). The device belongs to Group 8 (medical devices).
Applies in the countries of the European Union, as well as in Norway and Switzerland.

6 Delivery checklist

- 4 Vibration absorber
- 1 Double open-end wrench 17/19 mm
- 1 Hexagonal pin-type spanner
- 1 Lubricating grease for trunnions
- 1 Operating instructions
- 3 Program-data sheet
- 1 Information on how to remove the packing
- 1 Installation information sheet

The rotor(s) and associated accessories are included in the delivery in the quantity ordered.

7 Unpacking the centrifuge and installation

Unpack and install the centrifuge in accordance with the instruction leaflets AH5005XX and AH5005-01XX provided.

8 Initial operation

- According to the laboratory instrument standards EN / IEC 61010-2-020 an emergency switch to separate power supply in the event of a failure must be installed in the building electrical system. This switch has to be placed remote from the centrifuge, preferred outside of the room in which the centrifuge is installed or near by the exit of this room.
- According to the laboratory instrument standards IEC 61010-1 a switch to separate the power supply agreeable to IEC 60947-5-3 must be installed in the building electrical system. This switch has to be placed close to the centrifuge, must be easily reachable for the user and has to be characterized as a separator for this unit.
- Types 5005-90, 5005-08, 5005-88 and 5005-98 are permanently connected devices.
- Types 5005-80, 5005-90, 5005-08, 5005-88 and 5005-98 may only be installed by authorized skilled personnel.

Please note the attached mounting and installation instructions AH5005-02EN.

- **Position the centrifuge in a stable and level manner in a suitable place. During set-up, the required safety margin of 300 mm around the centrifuge is to be kept according to EN / IEC 61010-2-020.**

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

- Do not place any object in front of the ventiduct. Keep a ventilation area of 300 mm around the ventiduct.
- Centrifuge with PE connector:
If required, connect the PE connector at the back of the device with an additional medical potential equalization system.
- Check whether the mains voltage tallies with the statement on the type plate.
- For types 5005, 5005-50, 5005-80 and SK 48.06-30:
Connect the centrifuge with the connection cable to a standard mains socket. For connection ratings refer to Chapter "Technical specification".
- Switch on the mains switch in the building installation.
- Switch on the switch on the front panel. Switch position "I".
The following display follows:
1. Type of centrifuge, 2. The last rotor code detected by the rotor identification and the maximum speed of the rotor (n-max-Rotor), 3. The programme version, 4. **OPEN** / **OEFFNEN**.
- Open the lid.
The centrifuge data of the last used program or program 1 will be displayed.

9 Opening and closing the lid

9.1 Opening the lid


The lid can only be opened when the centrifuge is switched on and the rotor is at rest. If it cannot be opened under these circumstances, see the section on "Emergency release".

- Press the button  on the front panel. The lid unlocks via the motor, the illuminated pushbutton  extinguishes and in the display the symbol  changes to the symbol .


9.2 Closing the lid

Do not put your fingers between lid and housing.
Do not bang the lid shut.

- Place the lid and slightly press it down at the tubular handle. The locking action is effected by motor. The button  lights up and the symbol  in the display turns to the symbol .

10 Installation and removal of the rotor


- Clean the motor shaft (C) and the rotor drilling (A), and lightly grease the motor shaft afterwards. Dirt particles between the motor shaft and the rotor hinder a perfect seating of the rotor and cause an irregular operation.
- Place the rotor vertically on the motor shaft. The motor shaft dog (D) has to fit in the rotor slot (B). The alignment of the groove is labelled on the rotor.
- Tighten the rotor tension nut with the supplied wrench by turning in a clockwise direction.
- Check the rotor for firm seating.
- Loosening the rotor: Loosen the tension nut by turning in a counter clockwise direction, and turning until the working point for lifting. After passing the working point for lifting the rotor is loosened from the motor shaft cone. Turn the tension nut until the rotor is able to be lifted from the motor shaft.

11 Loading the rotor


Standard centrifuge containers of glass will not stand RCF values exceeding 4000 (DIN 58970, pg. 2).

- Check the rotor for firm seating.
With swing-out rotors all rotor positions must be lined with **identical** hangers. Certain hangers are marked with the number of the rotor position. These hangers may only be used in the respective rotor position. Hangers that are marked with a set number (e.g. S001/4) may only be used in the set.
- The rotors and hangers may only be loaded symmetrically. The centrifuge containers have to be distributed evenly on all rotor positions. For authorised combinations see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".
In the case of angle rotors all possible rotor positions must be loaded, see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

Rotor is evenly loaded

Not permitted!
Rotor is not evenly loaded

- On certain hangers, the weight of the maximum load or the weight of the maximum load and the maximum weight of the completely loaded hanger is specified. This weight may not be exceeded. In case of exception, see chapter "Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³". The weight specified for the maximum loading includes the total weight of adapter, frame, centrifuging container and content.
- In containers with rubber inserts, the same number of rubber inserts must always be among the centrifuge containers.
- Always fill the centrifuge containers outside of the centrifuge.
- No liquid should be allowed to enter the centrifugal chamber during filling and swinging out of the hangers.
- The maximum filling quantity for the centrifuge containers specified by the manufacturer must not be exceeded.
- In order to maintain the weight differences within the centrifuge container as marginal as possible, a consistent fill level in the containers is to be heeded.

12 Control and display devices

See figure on page 2.

Fig. 2: Display and operating panel

12.1 Status symbols

Lid open.

Lid closed.

Rotation indicator. The indication appears while the centrifuge is running as long as the rotor is turning.

STOP

Centrifugation run stopped or finished. Indication after completion of the centrifugation run as long as the rotor is turning. The indicator flashes following an emergency stop.

**LOCK 1,
LOCK 2**

Keyswitch position (for centrifuges with keyswitch only).

**LOCK 4,
LOCK 5**

Program interlocking when serial communication provided (for centrifuges with serial communication only).

PC, PE

Serial communication (for centrifuges with serial communication only).

Operating errors or faults that occur are shown symbolically in the display (see Chapter "Faults").

12.2 Control knob

For setting the individual parameters.

Turning anticlockwise reduces the value. Turning clockwise increases the value.

12.3 Keys and set-up options

Running time, parameter t/min:sec

1. Parameter t/min: Adjustable from 1 - 999 min, in 1 min steps.
2. Parameter t/ :sec Adjustable from 1 - 59 s, in 1 second steps.
3. Continuous running "---:--"

Interrogation of the integral RCF, parameter ∫RCF.

1. Speed, parameter RPM

A numerical value from 50 RPM up to the maximum speed of the rotor (n-max-Rotor) can be set. Maximum speed of the rotor, see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories". Adjustable in steps of 10.

2. Interrogation of the maximum speed of the rotor, parameter n-max-Rotor.

1. Relative centrifugal force, parameter RCF/RZB

A numerical value can be set, which gives a speed between 50 RPM and the maximum speed of the rotor (n-max-Rotor). Adjustable in steps of 1.

2. Interrogation of the maximum RCF of the rotor, parameter RCF-max-Rotor.

Run-up parameters

1. Run-up steps, parameter \curvearrowright_{1-9} . Step 9 = shortest run-up time, ... Step 1 = longest run-up time.
2. Run-up time, parameter $\curvearrowright_{min:sec}$. The time range that can be set is dependent upon the set speed.

Run-down parameters

1. Braking steps, parameter \curvearrowleft_{0-9} . R = linear braking curve, B = similar to an exponential braking curve. Step R9, B9 = short run-down time, ... Step R1, B1 = long run-down time, Step R0 = unbraked run-down.
2. Run-down time, parameter $\curvearrowleft_{min:sec}$. The time range that can be set is dependent upon the set speed.
3. Brake de-energisation speed, parameter $n^{(+)}/RPM$
After reaching this speed, unbraked run-down occurs.

1. Temperature setpoint, parameter T/°C (for centrifuges with cooling only). Adjustable from -20°C to +40°C, in 1°C steps (heating/cooling option settable between -20°C and +60°C / +90°C). The lowest temperature that can be achieved is dependent upon the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories").

2. Centrifuging radius, parameter r/mm. Entry in mm. Centrifuging radius, see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

1. Start centrifugation run. The rotation indicator  appears.
2. Acceptance of changes during the centrifugation run.

Stop centrifugation run.

The rotor decelerates with the preset run-down parameters. Pressing the key twice will initiate an emergency stop.

Select program location, parameter PROG-Nr.

Calling up of programs.

Saving of programs.

89 programs can be saved (program locations 1 to 89). Note: The program locations "----" and 90 to 99 are used as automatic buffers (see Chapter "Automatic buffer"). No programs can be saved in these program locations.

13 Entering centrifugation parameters

It is only possible to enter a parameter when the input field is displayed inversely (dark background). An inverse input field will go off automatically after 10 seconds.

13.1 Running time

13.1.1 Time selection

- Press the key repeatedly until the parameter t/min: or t/:sec is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

13.1.2 Continuous running

- Select the parameters t/min: and t/:sec one after the other (see Chapter "Time selection"), and set both to zero using the control knob . "---:--" will be displayed in the input field.

13.2 Speed (RPM)

- Press the key repeatedly until the parameter RPM is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

13.3 Relative centrifugal force (RCF/RZB)

- Press the key repeatedly until the parameter RCF/RZB is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

The centrifuging radius indicator will flash while the adjustment is being made.

13.4 Run-up and run-down parameters

The run-up and run-down parameters that have been set are displayed.

x: 1-9 = run-up step, t = run-up time

y: R1-R9, B1-B9 = braking step, R0 = unbraked run-down, t = run-down time, n^(*) = brake de-energisation speed

13.4.1 Run-up step

- Press the key repeatedly until the parameter ₁₋₉ or is displayed. The input field will be displayed inversely.
- Set the desired step by means of the control knob .

13.4.2 Run-up time

- Press the key repeatedly until the parameter _{min:sec} is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

If a run-up time is set that is longer than the running time, the centrifugation run will end before the set speed is reached.

13.4.3 Braking step

- Press the key repeatedly until the parameter ₀₋₉ or is displayed. The input field will be displayed inversely.
- Set the desired step by means of the control knob .

B-braking steps can only be set for special rotors.

13.4.4 Run-down time

If a brake de-energisation speed has been set, it is not possible to set a run-down time.

- Press the key repeatedly until the parameter _{min:sec} is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

13.4.5 Brake de-energisation speed

- Press the  key repeatedly until the parameter $n^{(4)}/RPM$ is displayed. The input field will be displayed inversely.
- Set the desired value by means of the control knob .

13.5 Radius/Temperature

13.5.1 Radius

- Press the  key repeatedly until the parameter r/mm is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

The value of the RCF/RZB adapts automatically when the radius is changed; this is shown by the indication flashing.

13.5.2 Temperature (for centrifuges with cooling only)

- Press the  key repeatedly until the parameter $T/^{\circ}C$ is displayed and the input field is shown inversely.
- Set the desired value by means of the control knob .

13.6 Automatic buffer

The buffer comprises the program locations "----" and 90 to 99. Every time a centrifugation run is started, modified centrifugation data are automatically stored in program location "----". The modified centrifugation data for the last 11 centrifugation runs are stored in the buffer and can be called up (see Chapter "Calling up programs").

14 Programming

14.1 Entering / modifying programs

- Set the desired parameters (see Chapter "Entering centrifugation parameters").
- Press the  key in order to select the parameter PROG-Nr. The input field will be displayed inversely.
- Set the desired program location by means of the control knob . If the program location indicator flashes, then this program location has already had centrifugation data assigned to it. In this case, set a free program location or overwrite the centrifugation data by continuing.
- Press the  key once to save the settings in the desired program location. Press the  key twice to overwrite centrifugation data already stored.

14.2 Calling up programs

- Press the  key in order to select the parameter PROG-Nr. The input field will be displayed inversely.
- Set the desired program location by means of the control knob .
- Press the  key. The centrifugation data for the selected program location will be displayed.

15 Centrifugation

When the centrifuge is running, according to EN / IEC 61010-2-020, no persons, dangerous substances or objects may be within the safety margin of 300 mm around the centrifuge.

If the permissible weight difference within the rotor loading has been exceeded, the drive shuts down during the start-up, and **IMBALANCE / UNWUCHT** is displayed.

If the speed in the selected program is higher than the maximum speed of the rotor ($n_{max-Rotor}$), it will not be possible to start a centrifugation run. The message **N > ROTOR MAX 96** will be displayed (see Chapter "Faults").

A centrifugation run can be stopped at any time by pushing the key .

It is possible to switch over between RPM and RCF display at any time using the  and  keys. When working with the RCF display, it is necessary to enter the centrifuging radius.

If **OPEN OEFFNEN** is displayed, a further operation of the centrifuge is only possible after opening the lid once.

If **R xx n-max xxxxx** is displayed, then no centrifugation run has taken place as the rotor was changed beforehand, refer to Chapter "Rotor Identification".


- Switch on the mains supply switch. Switch position **I**.
- Load the rotor and close the centrifuge lid.

15.1 Centrifugation with time selection

- Set the time or call up a program with a pre-set time (see Chapter "Entering centrifugation parameters" or "Calling up programs").
- Press the **START** key. The rotation indicator  will appear as long as the rotor is turning.
- When the time has expired or if the centrifugation run is interrupted by pressing the **STOP** key, run-down will take place with the selected run-down parameters. When the rotor comes to a stop, the  symbol will flash in the display until the lid is opened.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature (for centrifuges with cooling only) and the remaining time will be displayed.

15.2 Continuous running


- Select the **---:--** symbol or call up a continuous running program (see Chapter "Entering centrifugation parameters" or "Calling up programs").
- Press the **START** key. The rotation indicator  will appear as long as the rotor is turning. The time count starts from 00:00.
- Press the **STOP** key to stop the centrifugation run. Run-down will take place with the selected run-down parameters. When the rotor comes to a stop, the  symbol will flash in the display until the lid is opened.

During the centrifugation run, the speed of the rotor or the RCF value resulting from this, the sample temperature (for centrifuges with cooling only) and the time run will be displayed.

16 Changing settings during the centrifugation run.

The running time, the speed, the relative centrifugal force (RCF/RZB), the run-up and run-down parameters as well as the temperature (for centrifuges with cooling only) can be changed while the centrifuge is running.

The parameters can only be changed one at a time and in succession.

- Select the desired parameter and change the value using the control knob  (see Chapter "Entering centrifugation parameters").
- Press the **START** key. The modified setting will be stored in program location "----" (see Chapter "Automatic buffer"). The original program will not be overwritten.

17 Integral RCF (\int RCF)

The integral RCF (\int RCF) is a measure for the sedimentation effect ($\int n^2 dt$). This numerical value is used to compare centrifugation runs. To interrogate the integral RCF, press and hold the **\int RCF** key.

18 Displaying the maximum speed of the rotor

- Press the **\square** key repeatedly until the parameter RPM is displayed and the input field is shown inversely.
- Press the **\square** key once more and hold it pressed; the maximum speed of the rotor (n-max-Rotor) will be displayed.

19 Displaying the maximum RCF of the rotor

- Press the **RCF** key repeatedly until the parameter RCF/RZB is displayed and the input field is shown inversely.
- Press the **RCF** key once more and hold it pressed; the maximum RCF of the rotor (RCF-max-Rotor) will be displayed.

20 Emergency stop

- Press the **STOP** key twice. The **STOP** symbol will flash in the display.

In the case of an emergency stop, run-down takes place with braking step R9 (shortest run-down time). If braking step R0 was selected, the run-down time will be longer than with braking step R9 for technical reasons.

21 Audible signal

The audible signal sounds in accordance with the following philosophy:

OFF	<ul style="list-style-type: none"> at 2 s intervals on the occurrence of a fault
ON1	<ul style="list-style-type: none"> at 2 s intervals on the occurrence of a fault at 30 s intervals on completion of a centrifugation run and the rotor coming to rest.
ON2	<ul style="list-style-type: none"> at 2 s intervals on the occurrence of a fault at 30 s intervals on completion of a centrifugation run and the rotor coming to rest. every time a key is pressed

The audible signal is silenced by opening the lid or pressing any key.

The signal can be set as follows when the rotor is at standstill.

- Open the lid.
- Press and hold the **[I]** key for 8 seconds.
After 8 seconds **SOUND / BELL XXX** will be displayed.
- Set the desired function (OFF, ON1, ON2) by means of the control knob **○**.
- Press the **[START]** key to save the setting.
***** ok ***** will be displayed briefly as a confirmation.

22 Interrogating the operating hours

It is only possible to interrogate the operating hours when the rotor is at standstill.

- Open the lid.
- Press and hold the **[I]** key for 8 seconds.
After 8 seconds **SOUND / BELL XXX** will be displayed.
- Press the **[I]** key once more.
The operating hours (**CONTROL:**) for the centrifuge will be displayed.
The operating hours display will go off automatically after 10 seconds.

23 Setting the date and time

It is only possible to set the date and time when the rotor is at standstill.

- Open the lid.
- Press and hold the **[I]** key for 8 seconds.
After 8 seconds **SOUND / BELL XXX** will be displayed.
- Press the **[I]** key twice more.
The date and time will be displayed (a: year, mon: month, d: day, h: hours, min: minutes).
- Press the **[]** key repeatedly until the input field for the desired parameter is shown inversely.
- Set the desired value by means of the control knob **○**.
- Press the **[START]** key to save the setting.
***** ok ***** will be displayed briefly as a confirmation.
- To exit the date and time set-up, press any key apart from the **[]**, **[I]** and **[START]** keys.

24 Centrifugation data displayed after switching on.

After switching on, the centrifugation data for program 1 or for the program last used will be displayed. This can be set as follows:

- Open the lid.
- Switch the mains supply switch off and on again. Switch position **I**.
- Press the **[STOP]** key when the first visual change appears in the display (inverse display).
PROGRAM 1, LAST PROGRAM will be displayed.
- Set the desired function by means of the control knob **○**.
- Press the **[START]** key to save the setting.
***** ok ***** will be displayed briefly as a confirmation.

25 Immediate display of the centrifugation data after switching on

- Switch on the mains supply switch. Switch position **I**.
- Press any key apart from the **[STOP]** key when the first visual change appears in the display (inverse display). The centrifugation data will be displayed immediately.

26 Key-operated switch (for centrifuges with key-operated switch only).

The following program interlocks can be set by the key-operated switch:

- Left key position: **LOCK 1** will be displayed.
 Programs can be called up but not changed.
- Right key position: **LOCK 2** will be displayed.
 No programs can be called up or changed.
- Middle key position: no status display.
 No program interlock. Programs can be called up and changed.

27 Program linking (for centrifuges with program linking only).

Several centrifugation runs can be linked together with the help of the program linking function.

27.1 Linking programs

Program linking is only possible with programs in which run-up and braking steps have been set.
 Before linking, the programs must be stored in the desired sequence either by entering the program or by calling up the program (see Chapter "Programming").
 The program locations must be consecutive (e.g. program locations 10+11+12).

1. Press the **PROG** key in order to select the parameter PROG-Nr. The input field will be displayed inversely.
2. Set the program location for the first program (XX+) by means of the control knob \odot .
3. Press the **RCL** key. The centrifugation data for the selected program location will be displayed.
4. Press the **PROG** key twice in order to select the parameter PR-PART. The input field will be displayed inversely.
5. Press the **STO** key twice. The program will be linked and the program number for the next program location (+XX+) will be displayed.
6. Press the **RCL** key twice. The centrifugation data for the selected program location will be displayed.
7. Press the **STO** key twice. The program will be linked and the program number for the next program location (+XX+) will be displayed.
8. Keep repeating steps 6 and 7 until all programs have been linked.
9. Press the **PROG** key to exit. The program number for the last program (+XX) will be displayed.


27.2 Changing program links

- Call up the desired program (see Chapter "Calling up programs"), change the desired parameters (see Chapter "Entering centrifugation parameters"), and save the modified centrifugation data to the same program location (see Chapter "Entering / modifying programs"). As a result of saving, the program linking will be removed.
- Re-link the programs (see Chapter "Linking programs").

27.3 Centrifugation run with program linking

- Press the **PROG** key twice in order to select the parameter PR-PART. The input field will be displayed inversely.
- Set the program location for the first program (XX+) by means of the control knob \odot .
- Press the **RCL** key. The centrifugation data for the selected program location will be displayed.
- Press the **START** key. The rotation indicator \odot will appear as long as the rotor is turning.

The run-up and braking steps for the program linking will be displayed:

- When the time in the last program has expired, run-down will take place with the braking step of the last program. If the centrifugation run is interrupted by pressing the **STOP** key, run-down will take place with the braking step of the program running at the time.

27.4 Deleting program linking

- Press the **PROG** key in order to display the parameter PROG-Nr. The input field will be displayed inversely.
- Set the program location for the first program (XX+) by means of the control knob \odot .
- Press the **RCL** key. The centrifugation data for the selected program location will be displayed.
- Press the **PROG** key twice in order to display the parameter PR-PART. The input field will be displayed inversely.
- Press the **STO** key twice.
- Press the **PROG** key.

28 Cooling (only in centrifuges with cooling)

The temperature set-point can be adjusted from -20°C to +40°C. The lowest obtainable temperature is dependent on the rotor (see Chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"). If the actual temperature and the set temperature deviate more than 5 °C, the temperature value in the display will start to flash on and off.

28.1 Standby-cooling

With rotor standstill and closed lid the centrifugal chamber is cooled to the pre-selected temperature. The temperature set-point is shown in the display.

28.2 Pre-cooling the rotor

For rapid pre-cooling of the unloaded rotor and the accessories, it is recommended that the centrifuge is run with the continuous running settings and a speed of ca. 20% of the maximum rotor speed.

29 Heating (only on centrifuges with heating/cooling option)

During operation of the centrifuge the centrifuge chamber can be heated as required to the preselected temperature. The heating is switched off when the rotor is stationary.

Risk of burning! The surface temperature of the heating element in the centrifugal chamber of the centrifuge can reach 500°C / 932°F. Do not touch the heating element.

30 Relative centrifugal force (RCF)

The relative centrifugal force (RCF) is given as a multiple of the acceleration of gravity (g). It is a unit-free value and serves to compare the separation and sedimentation performance.

These values are calculated using the formula below:

$$\text{RCF} = \left(\frac{\text{RPM}}{1000} \right)^2 \times r \times 1,118 \quad \Rightarrow \quad \text{RPM} = \sqrt{\frac{\text{RCF}}{r \times 1,118}} \times 1000$$

RCF = relative centrifugal force

RPM = rotational speed (revolutions per minute)

r = centrifugal radius in mm = distance from the centre of the turning axis to the bottom of the centrifuge.

For more on the centrifugal radius see the chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories".

The relative centrifugal force (RCF) stands in relation to the revolutions per minute and the centrifugal radius.

31 Centrifugation of materials or mixtures of materials with a density higher than 1.2 kg/dm³

When centrifuging with maximum revolutions per minute the density of the materials or the material mixtures may not exceed 1.2 kg/dm³.

The speed must be reduced for materials or mixtures of materials with a higher density.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{\text{red}}) = \sqrt{\frac{1.2}{\text{Greater density [kg/dm}^3\text{]}}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, density 1.6 kg/dm³

$$n_{\text{red}} = \sqrt{\frac{1.2 \text{ kg/dm}^3}{1.6 \text{ kg/dm}^3}} \times 4000 \text{ RPM} = 3464 \text{ RPM}$$

In the exceptional case that the maximum loading indicated on the hanger is exceeded, the speed must also be reduced.

The permissible speed can be calculated using the following formula:

$$\text{Reduced speed } (n_{\text{red}}) = \sqrt{\frac{\text{maximum load [g]}}{\text{actual load [g]}}} \times \text{maximum speed [RPM]}$$

e.g.: maximum speed RPM 4000, maximum load 300 g, actual load 350 g

$$n_{\text{red}} = \sqrt{\frac{300 \text{ g}}{350 \text{ g}}} \times 4000 \text{ RPM} = 3703 \text{ RPM}$$

If in doubt you should obtain clarification from the manufacturer.

32 Rotor Identification

After every start of a centrifugation run the rotor utilised is identified.

After a change of rotor the drive switches off and the rotor code (R xx) as well as the maximum rotational speed (n-max=xxxxx) of the rotor are displayed.

 A further operation of the centrifuge is only possible after a single opening of the lid.
If, following a rotor change, the maximum speed of the rotor is less than the set speed, the speed is limited to the maximum speed of the rotor.

33 Emergency release

During a power failure the lid cannot be unlocked by motor. An emergency release has to be executed by hand.

 For emergency release disconnect the centrifuge from the mains. For permanently connected devices, switch off the mains switch to disconnect the mains supply to the device in the building installation.
Open the lid only during rotor standstill.

See figure on page 2.

- Switch off the switch on the front panel (switch position "0").
- Look through the window in the lid to be sure that the rotor has come to a standstill.
- Insert the hexagon socket head wrench horizontally into the drilling (Fig. 1, A) and carefully rotate by half a turn in clockwise direction until the lid can be opened.
- Pull the hexagon socket head wrench out of the drilling again.
- After turning the centrifuge on again, press the button  so that the motor-driven lid locking once again assumes the normal position (opened).

34 Maintenance and servicing

Pull the mains plug before cleaning. For permanently connected devices, switch off the mains switch to disconnect the mains supply to the device in the building installation.

Before any other cleaning or decontamination process other than that recommended by the manufacturer is applied, the user has to check with the manufacturer that the planned process does not damage the device.

- Centrifuges, rotors and accessories must not be cleaned in rinsing machines.
- They may only be cleaned by hand and disinfected with liquids.
- The water temperature must be between 20 – 25°C.
- Only detergents/disinfectants may be used which:
 - have a pH between 5 - 8
 - do not contain caustic alkalis, peroxides, chlorine compounds, acids and alkaline solutions
- In order to prevent appearances of corrosion through cleaning agents or disinfectants, the application guide from the manufacturer of the cleaning agent or disinfectant are absolutely to be heeded.

34.1 Centrifuge (housing, lid and centrifuging chamber)

34.1.1 Surface cleaning and care

- Clean the centrifuge housing and the centrifuging chamber regularly, using soap or a mild detergent and a damp cloth if required. For one thing, this services purposes of hygiene, and it also prevents corrosion through adhering impurities.
- Ingredients of suitable detergents:
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove the detergent residue by wiping with a damp cloth.
- The surfaces must be dried immediately after cleaning.
- In the event of condensation water formation, dry the centrifugal chamber by wiping out with an absorbent cloth.
- Lightly rub the rubber seal of the centrifuge chamber with talcum powder or a rubber care product after each cleaning.
- The centrifuging chamber is to be checked for damage once a year.

If damage is found which is relevant to safety, the centrifuge may no longer be put into operation. In this case, notify Customer Service.

34.1.2 Surface disinfection

- If infectious materials penetrates into the centrifugal chamber this is to be disinfected immediately.
- Ingredients of suitable disinfectants:
ethanol, n-propanol, isopropyl alcohol, glutardialdehyde, quaternary ammonium compounds.
- After using disinfectants, remove the disinfectant residue by wiping with a damp cloth.
- The surfaces must be dried immediately after disinfecting.

34.1.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for removing radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove the agent residue by wiping with a damp cloth.
- The surfaces must be dried directly after removing the radioactive contaminants.

34.2 Rotors and Attachments

34.2.1 Cleaning and care

- In order to prevent corrosion and material changes, rotors and accessories must be cleaned regularly with soap or a mild detergent and a damp cloth. Cleaning is recommended at least once a week. Contaminants must be removed immediately.
- Ingredients of suitable detergents:
soap, anionic tensides, non-ionic tensides.
- After using detergents, remove detergent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after cleaning.
- Angle rotors, container and hanger made of aluminium are to be lightly greased after drying using acid-free grease, e.g. vaseline.
- In the case of biosafety systems (for further details of available biosafety systems see chapter "Anhang/Appendix, Rotoren und Zubehör/Rotors and accessories"), the packing rings must be checked and cleaned regularly (weekly). The sealing ring is to be replaced immediately upon indication of crack formation, embrittlement or abrasive wear. In order to prevent the packing ring from twisting when opening and closing the cover, the packing ring must be lightly rubbed with talcum powder or a rubber care product.
- In order to prevent corrosion as a result of moisture between the rotor and the motor shaft, the rotor should be disassembled and cleaned at least once a month, and the motor shaft should be lightly greased.
- The rotors and the accessory parts have to be checked on a monthly basis for wear and corrosion. For swing-out rotors, it is important to check the area of the lifting lugs, for hangers, the grooves and the base should be checked for cracks.

Rotors and attachments may no longer be utilised upon indication of wear and tear or corrosion.

- Check the firm seating of the rotor on a weekly basis.

34.2.2 Disinfection

- If infectious material should get on the rotors or accessories, they must be appropriately disinfected.
- Ingredients of suitable disinfectants:
glutaraldehyde, propanol, ethyl hexanol, anionic tensides, corrosion inhibitors.
- After using disinfectants, remove disinfectant residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after disinfection.

34.2.3 Removal of radioactive contaminants

- The agent must be specifically labelled as being an agent for the removal of radioactive contaminants.
- Ingredients of suitable agents for removing radioactive contaminants:
anionic tensides, non-ionic tensides, polyhydrated ethanol.
- After removing the radioactive contaminants, remove agent residue by rinsing with water (only outside of the centrifuge) or wipe off with a damp cloth.
- The rotors and accessories must be dried directly after removing the radioactive contaminants.

34.2.4 Trunnions

With swing-out rotors the trunnions must be regularly lubricated (Hettich Lubricating Grease No. 4051) in order to ensure consistent swinging out of the hangers.

34.2.5 Rotors and accessories with limited service lives

The use of certain rotors, hangers and accessory parts is limited by time.

These are marked with the maximum permitted number of operating cycles or with an expiration date and the maximum permitted number of operating cycles or just with the expiration date; e.g.:

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For safety reasons, rotors, hangers and accessory parts may no longer be used if either the indicated maximum number of operating cycles or the indicated expiration date has been reached.

34.3 Autoclaving

Swing-out rotors, angle rotors made of aluminium, suspension made of metal, lids with bio-containment as well as stands and reductions can be autoclaved at 121° C / 250°F (20 mins).

Otherwise you must ask the manufacturer.

No statement can be made about the degree of sterility.

The lids of the rotors and containers must be removed prior to autoclaving.

Autoclaving accelerates the ageing process of plastics. In addition, autoclaving may discolour plastics.

We recommend that the packing rings of the bio-safety system be replaced after autoclaving.

34.4 Centrifuge containers

- With leakiness or after the breakage of centrifuging containers broken container parts and leaked centrifugation material are to be completely removed.
- The rubber inserts as well as the plastic sleeves of the rotors are to be replaced after a glass breakage.

Remaining glass splinters cause further glass breakage!

- If this concerns infectious material, a disinfection process is to be executed immediately.

35 Faults

If the fault cannot be eliminated with the help of the fault table, please inform Customer Service.

Please specify the type of centrifuge and the serial number. Both numbers can be found on the name plate of the centrifuge.

 Perform a MAINS RESET:

- Switch off the mains switch (switch position "0").
- Wait at least 10 seconds and then switch on the mains switch again (switch position "1").

Message / fault		Cause	Remedy
No display	---	No voltage. Overvoltage protection tripped out (only for types 5005-08).	<ul style="list-style-type: none"> - Check supply voltage. - Mains switch ON.
TACHO - ERROR	01	Faulty speedometer.	<ul style="list-style-type: none"> - Open the cover. - Switch off the mains switch (switch position "0"). - Wait at least 10 seconds. - Turn the rotor vigorously by hand. - Switch on the mains switch again (switch position "1"). The rotor must turn during switch-on.
	02	No rotor installed. Defective motor, frequency converter or drive.	
IMBALANCE / UNWUCHT	---	The rotor is unevenly loaded.	<ul style="list-style-type: none"> - Open lid. - Check the loading of the rotor, see chapter "Loading the rotor". - Repeat the centrifugation run.
CONTROL - ERROR	04, 06 - 09	Error in lid locking or lid closure.	<ul style="list-style-type: none"> - Perform a MAINS RESET.
N > MAX	05	Rotation too fast	
N < MIN	13	Rotation too slow	
ROTORCODE	10	Incorrect rotor coding	
MAINS INTERRUPT	---	Power failure during the centrifugation run. (The centrifugation run was not finished.)	<ul style="list-style-type: none"> - Open lid. - Push START button. - Repeat the centrifugation run if necessary.
VERSIONS-ERROR	12	Mismatch between electronic components	<ul style="list-style-type: none"> - Perform a MAINS RESET.
SER I/O - ERROR	30 - 38	Error / defective interface	
° C * - ERROR	50 - 56, 58	Error / defective cooling	
LOCK - ERROR	57	Error / defective program locking	
FU / CCI - ERROR	60 - 83	Error / defective motor control	
CONTROL-ERROR	26, 90 - 95, 97 - 99	Error / defective control unit	
N > ROTOR-MAX	96	Speed in the selected program greater than the maximum speed of the rotor (n-max-Rotor).	<ul style="list-style-type: none"> - Check the set speed. - Reduce the set speed

36 Acceptance of the centrifuges for repair

If the centrifuge is returned to the manufacturer for repair, it must be decontaminated and cleaned to protect persons, environment and material.

We reserve the right to accept contaminated centrifuges.

Costs incurred for cleaning and disinfection are to be charged to the customer.

We ask for your understanding in this matter.

37 Disposal

Before disposal, the device must be decontaminated and cleaned to protect people, the environment and property.

When you are disposing of the device, the respective statutory rules must be observed.

Pursuant to guideline 2002/96/EC (WEEE), all devices supplied after August 13, 2005 may not be disposed as part of domestic waste. The device belongs to group 8 (medical devices) and is categorized in the business-to-business field.

The icon of the crossed-out trash can shows that the device may not be disposed as part of domestic waste.

The waste disposal guidelines of the individual EC countries might vary. If necessary, contact your supplier.

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1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature

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



























1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
 2) Centrifuge containers of glass will not stand RCF values exceeding 4000












4176		4522							
Ausschwingrotor 6-fach / Swing out rotor 6-times 									
		4226	4227	4232	4233		4238	4241	4245
									
	Microtainer	0528	0509	4234	0551	5127		0513	
									
Kapazität / capacity	ml	0,8	100	15	750	600	250	25	50
Maße / dimensions \varnothing x L	mm	8 x 45	58 x 161	17 x 120	93 x 137	93 x 134	62 x 122	25 x 90	29 x 115
Anzahl p. Rotor / number p. rotor		468	6	138	6	6	6	48	48
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	²⁾	5818	5818	5999	5999	5999	5818	5818	5999
Radius / radius	mm	257	257	265	265	265	257	257	265
 9 (97%)	sec							118	
 9	sec							180	
Temperatur / temperature	°C ¹⁾							14	






4176		4522				4514	4524			
Ausschwingrotor 6-fach / Swing out rotor 6-times 										
		4249	4258	6322	SK 09.00	Kunststoff-Einsatz / Synth. mat. ins.	4529-O ¹⁰⁾	4529-M ¹⁰⁾	4529-U ¹⁰⁾	
										
		0512	Corning							
										
Kapazität / capacity	ml	50	750	250	14	100	500			
Maße / dimensions \varnothing x L	mm	29 x 115	97 x 152	60 x 172	16,5 x 106	---	---			
Anzahl p. Rotor / number p. rotor		36	6	6	138	18	12			
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500			
RZB / RCF	²⁾	5909	5999	5818	5954	6181	6498			
Radius / radius	mm	261	265	257	263	273	287			
 9 (97%)	sec							118		
 9	sec							180		
Temperatur / temperature	°C ¹⁾	14						16		

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







	Position of the holding pins	Distance of the holding pins from the bottom of the inserts (inside)
4529-O	upper section	199.5 mm
4529-M	middle section	182.0 mm
4529-U	lower section	164.5 mm







4176	4524			4546	4591	4595		
Ausschwingrotor 6-fach / Swing out rotor 6-times  								
	Baxter/ 4525	Biopack/ 4525	Fenwal 4525 6)	4592-H	4559-H	4592-H		4596
								
	2 x 4526	2 x 4527	2 x 4528				0550	4597
							 ⁸⁾	
Kapazität / capacity ml	500	500	500	500	450	500	2000	100
Maße / dimensions Ø x L mm	----	----	----	----	----	----	180x150x100	----
Anzahl p. Rotor / number p. rotor	12	12	12	12	12	12	6	12
Drehzahl / speed RPM	4500	4500	4500	4500	4500	4500	3500	1000
RZB / RCF	6498	6498	6498	6498	6271	6498	3848	312
Radius / radius mm	287	287	287	287	277	287	281	279
 9 (97%) sec							118	
 9 sec							180	
Temperatur / temperature °C ¹⁾							16	
							-10	

4176	4579	4579 + 4255 / 4255-P 4)						
Ausschwingrotor 6-fach / Swing out rotor 6-times  								
			4449	4430				
								
		4255 / 4255-P 4)	0512	4239	Corning		Nalgene®	Nunc®
		 ^{5) 8)}	 ^{5) 8)}	 ⁵⁾	 ⁵⁾	 ⁵⁾	 ⁵⁾	
Kapazität / capacity ml	1000	750	1000	500	250	175	200	
Maße / dimensions Ø x L mm	98 x 138	97 x 152	98 x 175	96 x 147	60 x 172	61,5 x 144	60 x 130	
Anzahl p. Rotor / number p. rotor	6	6	6	6	6	6	6	
Drehzahl / speed RPM	4500	4500	4500	4500	4500	4500	4500	
RZB / RCF	6294	6294	6294	6294	6294	6294	6294	
Radius / radius mm	278	278	278	278	278	278	278	
 9 (97%) sec							118	
 9 sec							180	
Temperatur / temperature °C ¹⁾							4	

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





- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 5) 4255 cannot be closed with the lid
- 6) as soon as Macopharma, Terumo, Kawasumi
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.







4176		4579 + 4255 / 4255-P 4)								
Ausschwingrotor 6-fach / Swing out rotor 6-times  ↙ 90°										
		4432		4433				4434		
										
		2078	0536		0578			0500	2079	0518
										
Kapazität / capacity	ml	1,5	2,0	5	7	2,7 - 3	4,5 - 5	9	10	15
Maße / dimensions	∅ x L mm	11 x 38	12 x 75	12 x 75	12 x 100	11 x 66	11 x 92	14 x 100	17 x 70	17 x 100
Anzahl p. Rotor / number p. rotor		252	180	180	180	180	180	114	114	114
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	²⁾	5796	5750	5750	5750	5750	5750	5886	5886	5886
Radius / radius	mm	256	254	254	254	254	254	260	260	260
 9 (97%)	sec									118
 9	sec									180
Temperatur / temperature	°C ¹⁾									4

4176		4579 + 4255 / 4255-P 4)								
Ausschwingrotor 6-fach / Swing out rotor 6-times  ↙ 90°										
		4434								
										
		0507								
										
Kapazität / capacity	ml	15	10	8	4 - 5,5	9 - 10	4 - 7	8,5 - 10	10	10
Maße / dimensions	∅ x L mm	17 x 100	16 x 80	16 x 81	15 x 75	16 x 92	16 x 75	16 x 100	15 x 102	15 x 102
Anzahl p. Rotor / number p. rotor		114	114	114	114	114	114	114	114	114
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	²⁾	5886	5886	5886	5886	5886	5886	5886	5886	5886
Radius / radius	mm	260	260	260	260	260	260	260	260	260
 9 (97%)	sec									118
 9	sec									180
Temperatur / temperature	°C ¹⁾									4

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





- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 5) 4255 cannot be closed with the lid








4176		4579 + 4255 / 4255-P 4)									
Ausschwingrotor 6-fach / Swing out rotor 6-times  ∠ 90°											
										4435	
						0509				0519	
											
Kapazität / capacity	ml	2,6 – 2,9	4,9	1,6 - 5	4 – 7	15	25	30	25		
Maße / dimensions	∅ x L	mm	13 x 65	13 x 90	13 x 75	13 x 100	17 x 120	25 x 90	25 x 110	24 x 100	
Anzahl p. Rotor / number p. rotor		126	126	126	126	72	42	42	42		
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500		
RZB / RCF	2)	5750	5750	5750	5750	6022	5818	5818	5615		
Radius / radius	mm	254	254	254	254	266	257	257	248		
	9 (97%)	sec								118	
	9	sec								180	
Temperatur / temperature	°C 1)									4	

4176		4579 + 4255 / 4255-P 4)															
Ausschwingrotor 6-fach / Swing out rotor 6-times  ∠ 90°																	
																4439	
				0521		Falcon®		Falcon®		0513		0526		5127		0528	
																	
Kapazität / capacity	ml	50	225	175	50	100	250	450	100								
Maße / dimensions	∅ x L	mm	34 x 100	61 x 137	61 x 118	29 x 115	44 x 100	62 x 122	70 x 159	58 x 161							
Anzahl p. Rotor / number p. rotor		24	6	6	30	12	6	6	6								
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500								
RZB / RCF	2)	5705	6294	6294	5999	5683	6113	6294	6294								
Radius / radius	mm	252	278	278	265	251	270	278	278								
	9 (97%)	sec								118							
	9	sec								180							
Temperatur / temperature	°C 1)									4							

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



- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature (only with cooling centrifuges)
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 3) tube will not stand RCF values exceeding 700
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 5) 4255 cannot be closed with the lid
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.



4176		4579 + 4255 / 4255-P 4)						
Ausschwingrotor 6-fach / Swing out rotor 6-times  								
		4466	SK 39.93	SK 61.98				
								
		0551	0530					
	 8)							
Kapazität / capacity	ml	600	250	50				
Maße / dimensions $\varnothing \times L$	mm	93 x 134	65 x 115	29 x 115				
Anzahl p. Rotor / number p. rotor		6	6	30				
Drehzahl / speed	RPM	4500	4500	4500				
RZB / RCF	2)	6294	5841	5796				
Radius / radius	mm	278	258	256				
 9 (97%)	sec		118					
 9	sec		180					
Temperatur / temperature	°C 1)		4					



4176		4547 + 5621								
Ausschwingrotor 6-fach / Swing out rotor 6-times  										
		mit Bioabdichtung / with bio-containment 9)								
				4449	4430			4432		
										
		0512	4239	Corning	Corning	Nalgene	Nunc	2078	0536	
		 8)	 8)							
Kapazität / capacity	ml	750	1000	500	250	175	200	1,5	2,0	
Maße / dimensions $\varnothing \times L$	mm	97 x 152	96 x 176	96 x 147	60 x 162	61,5 x 144,3	60 x 130	11 x 38	11 x 38	
Anzahl p. Rotor / number p. rotor		6	6	6	6	6	6	252	252	
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500	
RZB / RCF	2)	6384	6384	6384	6384	6384	6348	5886	5886	
Radius / radius	mm	282	282	282	282	282	282	260	260	
 9 (97%)	sec					112				
 9	sec					195				
Temperatur / temperature	°C 1)					9				

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 4) 4255-P: special surface treatment for highest hygienic requirements
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".











<p>4176</p> <p>Ausschwingrotor 6-fach / Swing out rotor 6-times</p>  <p>∠ 90°</p>	<p>4547 + 5621</p> 									
	<p>mit Bioabdichtung / with bio-containment 9)</p>									
	<p>4433</p> 				<p>4434</p> 					
	0553	0578			0500	2079	0507			
										
Kapazität / capacity	ml	5	7	2,7 - 3	4,5 - 5	9	10	15	10	
Maße / dimensions	∅ x L	mm	12 x 75	12 x 100	11 x 66	11 x 92	14 x 100	17 x 70	17 x 100	16 x 80
Anzahl p. Rotor / number p. rotor		180	180	180	180	114	114	114	114	
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500	
RZB / RCF	²⁾	5841	5841	5841	5841	5977	5977	5977	5977	
Radius / radius	mm	258	258	258	258	264	264	264	264	
 9 (97%)	sec	112								
 9	sec	195								
Temperatur / temperature	°C ¹⁾	9								













<p>4176</p> <p>Ausschwingrotor 6-fach / Swing out rotor 6-times</p>  <p>∠ 90°</p>	<p>4547 + 5621</p> 								
	<p>mit Bioabdichtung / with bio-containment 9)</p>								
	<p>4434</p> 								

									
Kapazität / capacity	ml	8	4 - 5,5	7,5 - 8,2	9 - 10	4 - 7	8	8,5 - 10	12
Maße / dimensions	∅ x L	mm	16 x 81	15 x 75	15 x 92	16 x 92	16 x 75	16 x 125	16 x 101
Anzahl p. Rotor / number p. rotor		114	114	114	114	114	114	114	114
Drehzahl / speed	RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	²⁾	5977	5977	5977	5977	5977	5977	5977	5977
Radius / radius	mm	264	264	264	264	264	264	264	264
 9 (97%)	sec	112							
 9	sec	195							
Temperatur / temperature	°C ¹⁾	9							

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





- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".





4176		4547 + 5621								
Ausschwingrotor 6-fach / Swing out rotor 6-times  		 mit Bioabdichtung / with bio-containment ⁹⁾								
		4434		4435				4437	4438 + 0726	
										
		0518				0509		0519		
										
Kapazität / capacity	ml	10	15	2,6 – 2,9	4,9	1,6 – 5	4 – 7	15	25	
Maße / dimensions	∅ x L	mm	15 x 102	17 x 100	13 x 65	13 x 90	13 x 75	13 x 100	17 x 120	24 x 100
Anzahl p. Rotor / number p. rotor			114	114	126	126	126	72	42	
Drehzahl / speed	RPM		4500	4500	4500	4500	4500	4500	4500	
RZB / RCF	²⁾		5977	5977	5841	5841	5841	6090	5728	
Radius / radius	mm		264	264	258	258	258	262	253	
 9 (97%)	sec								112	
 9	sec								195	
Temperatur / temperature	°C ¹⁾								9	

4176		4547 + 5621								
Ausschwingrotor 6-fach / Swing out rotor 6-times  		 mit Bioabdichtung / with bio-containment ⁹⁾								
		4438		4439		4440		4441	4442	4443
										
				0521		0513		0526	5127	
									⁸⁾	
Kapazität / capacity	ml	25	30	50	225	175	50	100	250	
Maße / dimensions	∅ x L	mm	25 x 90	25 x 110	34 x 100	61 x 137	61 x 118	29 x 115	44 x 100	62 x 122
Anzahl p. Rotor / number p. rotor			42	42	24	6	6	30	12	6
Drehzahl / speed	RPM		4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF	²⁾		5909	5909	5773	6384	6348	6090	5750	6203
Radius / radius	mm		261	261	255	282	282	269	254	274
 9 (97%)	sec								112	
 9	sec								195	
Temperatur / temperature	°C ¹⁾								9	

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 8) At temperatures above 40 °C and/or poor filling of the tubes, these can go out of shape.
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

4176	4547 + 5621				4523			
Ausschwingrotor 6-fach / Swing out rotor 6-times  								
	mit Bioabdichtung / with bio-containment 9)							
	4443	4466	SK 39.93	SK 61.98	4516			
								
	 8)	 8)			Baxter	Biopack	Fenwal	4-fach / 4-times
Kapazität / capacity ml	290	600	250	50	500	500	500	500
Maße / dimensions Ø x L mm	62 x 137	93 x 134	65 x 115	29 x 115	---	---	---	---
Anzahl p. Rotor / number p. rotor	6	6	6	50	6	6	6	6
Drehzahl / speed RPM	4500	4500	4500	4500	4500	4500	4500	4500
RZB / RCF 2)	6203	6384	6158	5886	6359	6359	6359	6359
Radius / radius mm	274	282	272	260	281	281	281	281
 9 (97%) sec					112			
 9 sec					195			
Temperatur / temperature °C 1)					9			

4176	SK 06.07 + 5629						
Ausschwingrotor 6-fach / Swing out rotor 6-times  							
	mit Bioabdichtung / with bio-containment 9)						
	4626						
							
	MTP	MTP	MS	CP	DWP	Microtest-platten / plate Terasaki	QP
							
Maße / dimensions Ø x L mm	86x128x15	86x128x17,5	86x128x46	86x128x22	86x128x44,5	59x84x11	86x128x 83
Anzahl p. Rotor / number p. rotor	36	30	6	24	6	12	6
Drehzahl / speed RPM	4500	4500	4500	4500	4500	4500	4500
RZB / RCF 2)	5139	5139	5139	5139	5139	5139	5139
Radius / radius mm	227	227	227	227	227	227	227
 9 (97%) sec					112		
 9 sec					195		
Temperatur / temperature °C 1)					9		

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

















CP Kulturplatte /
Culture plate

DWP Deep Well Platte /
Deep well plate

MS Micronic System /
Micronic system

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- 1) Lowest possible temperature during maximum speed, 1 h running time and 20°C ambient temperature
- 2) Centrifuge containers of glass will not stand RCF values exceeding 4000
- 7) only executable in centrifuges type SK 48.06-30
- 9) In conformity with DIN EN 61010, part 2 – 020. Observe the notes for bio safety systems in chapters "Notes on safety" and "Maintenance and servicing".

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All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.