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کد محصول	شرح
IAP-0500	AutoCAT 2 WAVE 1 (IABP) system includes: • FiberOptix sensor technology • WAVE algorithm—proprietary timing software • Aortic flow timing method • AutoPilot mode of operation
IAP-0400	AutoCAT 2 WAVE 1 (IABP) system includes: • AutoPilot operational mode with signal software

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- FiberOptix[®] capable:
 - AP signal transmitted at speed of light
- Proprietary WAVE[®] algorithm
- Proprietary Aortic Flow Timing Method
- AutoPilot[®] mode of operation
- Microprocessor-based system architecture
- Modular system consisting of display/control module and pneumatic drive unit
- Proprietary deflation timing management

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- AC requirements:
 - 90–264 VAC 47–63 Hz
- Typical power consumption: 245 watts
- Maximum power consumption: 420 watts
- Battery operating time:
 - 90 minutes minimum with full charge
 - 180 minutes with optional second battery
- Typical battery recharging time:
 - 80% in 4 hours from full discharge
 - Recharge to 80% indicated by yellow light

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- Control module with monitor:
 - 10" high (25.4 cm) x 13.75" wide (35 cm) x 2" deep (5 cm)
- Pneumatic drive unit:
 - 31.5" high (80 cm) x 13.5" wide (34.3 cm) x 21" deep (53.3 cm)

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- Control module:
 - 5 lbs (2.3 kg)
- Pneumatic unit for AutoCAT 2 WAVE:
 - 95.5 lbs (42.4 kg)
- Total weight for AutoCAT 2 WAVE:
 - 100.5 lbs (44.7 kg)
- Total weight for AERO[®] Series:
 - 91.5 lbs (40.7 kg)

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- Drive system: Stepper motor-driven bellows
- Drive gas: USP-grade helium
- Helium tank:
 - Disposable canister (500 psi) or refillable (2000 psi) cylinder—US Approval; (2900 psi) cylinder—European Approval
- Pumping volume:
 - 0.5cc to 50cc, adjustable in 0.5cc increments
- Counterpulsation rate: 40 to 200 pulsations/minute
- Assist ratio options

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- Thermoelectric system removes moisture continuously from pneumatic system without interrupting counterpulsation

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- AutoPilot:
 - Automatically selects ECG/AP signal, sources, trigger mode, and timing method as well as timing settings
 - Automatically changes settings to optimize assist
 - Proprietary software sets timing to correspond to individual patient needs
- Operator:
 - Allows user control of most pump functions

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- ECG (PATTERN, PEAK, AFIB):
 - Microprocessor-based R-waveform trigger detection algorithms
- Pacer (VPACE, APACE):
 - Low level (skin) ECG input
 - Pulse width 0.1 to 0.5 ms and pulse amplitude => +5 to +700 mV
 - Pulse width => 0.5 to 2 ms and pulse amplitude => +2 to +700 mV
 - High level (monitor) input
 - Pulse width 0.1 to 2 ms and pulse amplitude => 1 V
 - AV pacer detection is <250 msec between pacer pulses
- Arterial pressure (AP):
 - Microprocessor-based waveform trigger detection algorithm
- Internal:
 - Default to 80 bpm; adjustable 40 to 120 bpm
- Filtering:
 - Diathermy, 30 Hz low pass

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
ECG TRIGGER MODES:

- PATTERN: HR <130 bpm no arrhythmia
- PEAK: HR >130 bpm or arrhythmia detected and arrhythmia timing OFF
- AFIB: Any HR with arrhythmia detected
- VPACE: Single or dual pacer (<250 msec apart) and no QRS or AP waveform detected
- APACE: Single pacer with R-wave >100 msec later. Transition only

AP TRIGGER MODE:

- No ECG signal or noisy ECG signal

*Based upon deflation timing management.

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INFLATION TIMING METHODS:

- Aortic Flow:
 - Proprietary WAVE algorithm sets the timing intra-beat on average 1 ms of aortic valve closure
- Predictive: AP waveform analysis to set inflation
Weissler: ECG only, inflation timing based on systolic time intervals

DEFLATION TIMING METHODS:

- R-wave: Real-time deflation on R-wave
- Predictive: Deflation set to occur just prior to next systolic rise
- Weissler: ECG only, deflation timing based on diastolic intervals

MANUAL:

- User set inflation and deflation timing in Operator Mode

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- ECG:
 - Inflation, 20%–80% of R-R interval
 - Deflation, 30%–120% of R-R interval
- AP:
 - Inflation, 0–35% of peak systole-peak systole interval
 - Deflation, 35%–75% of peak systole-peak systole interval
- AFIB Trigger Mode:
 - Inflation 80 to 430 ms after R-wave trigger event
 - Deflation on R-wave

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- Type:
 - Color LCD flat screen
- Channels:
 - Three-channel multicolor waveforms
 - ECG: Green trace with white highlight on assisted portion
 - AP: Red trace calibrated for direct reading of AP, white highlight on assisted portions when in Operator Mode
 - Balloon pressure: Blue trace calibrated in mm Hg and displayed continuously
- Timing reference display:
 - Numerical timing settings in both operating modes as well as a bar graph displaying inflate/deflate events in Operator Mode
- Cursor:
 - Measurement of AP and balloon pressure waveforms

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- Patient hemodynamics: Heart rate, AP—systolic, augmented, diastolic, and mean arterial. When in 1:2 or lower assist ratio the assisted values are displayed in white and the unassisted values are displayed in yellow
- Displayed parameters: ECG source and gain state, alarm status with timer, ON Battery indication, operation mode selection, AP alarm parameter and limit, timing settings, helium tank level, arrhythmia detection, and timing status
- Operations status: Operational mode, trigger mode, helium tank gauge, alarm/battery charge status, balloon volume
- Diagnostic alarm/help messages: Preprogrammed troubleshooting prompts/help

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- Recorder:
 - Dual-channel dot matrix: Dot density 400 dots/inch, 25 mm/s
- Waveforms:
 - ECG, AP, or balloon pressure (one or two recorded)
- Alphanumeric:
 - Operational mode, trigger mode, ECG lead/source, AP source, AP alarm status, timing settings, assist ratio, balloon volume, timing method, arrhythmia status, alarm condition, date, time, patient hemodynamics

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- Freezes approximately 7 seconds of patient data on screen

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- ECG:
 - 5 lead skin cable (I, II, III, aVR, aVL, aVF and V)
 - High level monitor input (0 to 5 V)
- AP:
 - Fiber optic signal input from LightWave[®]AB Catheter (WAVE)
 - AP transducer (Spectramed or equivalent), 50 mV/V/cm Hg
 - High-level monitor input (1 V = 100 mm Hg)

Note: Additional system specifications are available from Arrow upon request. Specifications are subject to change without notice.

Caution: U.S. Federal law limits this device to sale by or on order of a physician.

U.S. Patent Nos. 6,258,035, 6,569,103, 6,887,206; additional patents pending.

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WECK HORIZON METAL LIGATION SYSTEM

□ Micro ■ Small ■ Small-Wide ■ Medium ■ Medium- Large ■ Large

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WECK® HORIZON™ METAL LIGATION SYSTEM

□ Micro ■ Small ■ Small-Wide ■ Medium ■ Medium-Large ■ Large

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شماره کاتالوگ	طول	سایز کلیپس	نوع فک
537061	6" (15 cm)	MicroClip™	Curved
537081	8" (20 cm)	MicroClip™	Curved
137061	6" (15 cm)	Small	Curved
137081	8" (20 cm)	Small	Curved
137085	8" (20 cm)	Small	Angled
137111	11" (28 cm)	Small	Curved
137062	6" (15 cm)	Small-Wide	Curved
137082	8" (20 cm)	Small-Wide	Curved
137086	8" (20 cm)	Small-Wide	Angled
137112	11" (28 cm)	Small-Wide	Curved
237061	6" (15 cm)	Medium	Curved
237081	8" (20 cm)	Medium	Curved
237085	8" (20 cm)	Medium	Angled
237111	11" (28 cm)	Medium	Curved
237115	11" (28 cm)	Medium	Angled
237117	11" (28 cm)	Medium	Right-Angled
237141	14" (35.5 cm)	Medium	Curved
237145	14" (35.5 cm)	Medium	Angled
337081	8" (20 cm)	Medium-Large	Curved
337085	8" (20 cm)	Medium-Large	Angled
337111	11" (28 cm)	Medium-Large	Curved
337115	11" (28 cm)	Medium-Large	Angled
337117	11" (28 cm)	Medium-Large	Right-Angled
337141	14" (35.5 cm)	Medium-Large	Curved
337145	14" (35.5 cm)	Medium-Large	Angled
437081	8" (20 cm)	Large	Curved
437085	8" (20 cm)	Large	Angled
437111	11" (28 cm)	Large	Curved
437115	11" (28 cm)	Large	Angled
437117	11" (28 cm)	Large	Right-Angled
437141	14" (35.5 cm)	Large	Curved
437145	14" (35.5 cm)	Large	Angled

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