


Abadis Med



Filters

www.abadis-med.com



Antibacterial filter



- Easy to use.
- It has 99/99% filtration capability.
- can filter microorganisms up to 0.3 microns.
- It is used in all types of portable and hospital suction.
- Prevents the spread of microorganisms to the environment.

Filter life span: When the paper inside the filter changes color, it should be replaced.

Flow stopper filter



- It is used in all types of portable and hospital suction.
- Used to prevent overflow of suctioned fluids.
- Used to prevent the spread of infectious fluids
- Prevents the release of toxic gases in electrosurgery and smoke.
- Prevents damage of the suction device and increases lifespan of suction device.

Filter life span: When liquids penetrate the filter surface or if suction pressure drops occur, the filter must be replaced.

Address: No. 104, Andarzgoo Blvd, Qeytariyeh area, Tehran, Iran

Tel: 021-22671845

Fax: 021-22219234