




شرکت مهندسی پارس نهند  
PARS NAHAND ENGINEERING Co.



## Parsonic 30s Ultrasonic Cleaner



|                                  |                                |
|----------------------------------|--------------------------------|
| Operating Voltage                | 220 V 50 Hz                    |
| Body, Tank, Basket & Drain Valve | Stainless Steel 304            |
| Tank Capacity & Dimension        | 28 Lit. 500x300x200h mm        |
| Body Dimension                   | 600x325x360h mm                |
| Ultrasonic Power & Frequency     | 400 Watt Max , 28 Khz ±5%      |
| Heater                           | 600 Watt 20-70 °C , Step 10 °C |
| Timer                            | 2-22 min , Step 4 min          |
| Accessory                        | Cable, Lid, Basket, Drain Hose |
| Guarantee                        | One Year                       |

ISO 9001:2015 & 13485:2016 Certificate

Iran Ministry of Health Certificate for Medical Application

تهران ، کیلومتر ۲۰ جاده آبعلی ، پارک فن آوری پردیس ، خیابان نوآوری پنجم  
پلاک ۵۱ کدپستی: ۱۶۵۴۱۲۰۸۷۷      تلفن / فاکس: ۵ - ۰۲۱ - ۷۶۲۵۰۱۸۱  
No. 51 , 5th Noavari Str. , Pardis TechPark , Tehran , Iran  
Zip Code: 1654120877      TEL / FAX: +98 21 76250181-5  
[www.pnec.co.ir](http://www.pnec.co.ir)      [info@pnec.co.ir](mailto:info@pnec.co.ir)

