

LAPAROSCOPY SYSTEM



LAPAROSCOPY SYSTEM

Endoset Flow110

- Max 40 L/min & 30 mmHg, with more sensitivity in Pediatric mode
- Display of gas consumption and bottle gas pressure
- intelligent Algorithm for fast & safe insufflation
- Built in pre gas-heating technology

Endoset Light110

- Lamp life >30000 h
- Light cable recognition
- Variable light intensity
- Homogeneous illumination and flicker free

Endoset ProLight110

- Lamp life >30000 h
- Broad band light spectrum with Variable light intensity
- Variant light cable support and light cable recognition
- Homogeneous illumination and flicker free

Endoset Video110

- Real Time Processing (WB, AE, AGC, Sharpness, 2X D-Zoom, ...)
- High resolution and 5 output (2*DVI, 2*SDI, 1*S-Video)
- Ergonomic and autoclavable camera head
- Image capture & video recorder (Optional)

Endoset View110

- For use in laparoscopic surgery (30cm, 10mm, 0 Degree)
- Autoclavable and ethylene oxide sterilization
- Fiber optic light transmission and Extra low dispersion glass

